

3. PCR Consumables


PCR 96 Well Plate And Sealing Film



Features:

- ✧ Made from Virgin PP to ensure the minimal loss of reaction solution
- ✧ Uniform thickness, better thermal transfer ensured
- ✧ Flat surface, thick and solid, not easy to deform
- ✧ Non-pyrogenic, Non- endotoxin, DNase and RNase free
- ✧ Fit with 100 μ L or 200 μ L heating block
- ✧ No skirt, semi-skirt and full-skirt


PCR 384 Well Plate



Features:

- ✧ Made from Virgin PP to ensure the minimal loss of reaction solution
- ✧ Flat surface, thick and solid, not easy to deform
- ✧ Elevated edge of the hole can better prevent cross contamination
- ✧ Compatible with sealing film, and fit with most common instruments
- ✧ No skirt and full-skirt

96 Deep Well Plate & Tip Comb



- ✧ Used for Nucleic acid extraction
- ✧ Non-pyrogenic, DNase and RNase free
- ✧ 96 well plate can be used for magnetic bead kit

Order Information (S series and J series with different delivery time and price)

S Series PCR Consumables

PCR 96-well Plate

Picture	Cat#	Description	Packaging (pcs/pk)
	26031001	0.1mL PCR 96 Well Plate, No-Skirt, Transparent, Cut Corner H12, without Printing, 10pcs*5bags/Carton	50
	26031003	0.1mL PCR 96 Well Plate, No-Skirt, Transparent, Black Printing, Cut Corner H12, 10pcs*5bags/Carton	50
	26031002	0.1mL PCR 96 Well Plate, No-Skirt, White, Cut Corner H12, without Printing, 10pcs*5bags/Carton	50
	26031004	0.1mL PCR 96 Well Plate, No-Skirt, White, Black Printing, Cut Corner H12, 10pcs*5bags/Carton	50
	26031005	0.1mL PCR 96 Well Plate, Half-Skirted, Transparent, Cut Corner H12, without Printing, 10pcs*5bags/Carton	50

	26031006	0.1mL PCR 96 Well Plate, Half-Skirted, White, Cut Corner H12, without Printing, 10pcs*5bags/Carton	50
	26031007	0.1mL PCR 96 Well Plate, Half-Skirted, Transparent, For ABI 7500 Fast, Cut Corner A1, without Printing, 10pcs*5bags/Carton	50
	26031014	0.2mL PCR 96 Well Plate, Half-Skirted, White, For ABI 7500, Cut Corner A12, without Printing, 10pcs*5bags/Carton	50
	26031009	0.2mL PCR 96 Well Plate, No-Skirt, Transparent, Cut Corner H1, without Printing, 10pcs*5bags/Carton	50
	26031012	0.2mL PCR 96 Well Plate, No-Skirt, Transparent, Black Printing, Cut Corner H1, 10pcs*5bags/Carton	50
	26031010	0.2mL PCR 96 Well Plate, No-Skirt, White, Cut Corner H1, without Printing, 10pcs*5bags/Carton	50
	26031013	0.2mL PCR 96 Well Plate, No-Skirt, White, Black Printing, Cut Corner H1, 10pcs*5bags/Carton	50
	26031008	0.2mL PCR 96 Well Plate, Half-Skirted, Transparent, Cut Corner A12, without Printing, 10pcs*5bags/Carton	50
	26031011	0.2mL PCR 96 Well Plate, Half-Skirted, Transparent, Black Printing, Cut Corner A12, 10pcs*5bags/Carton	50
	26041001	Pressure Sealing Film For PCR Plates, Transparent, Universal For 96 and 384 Well Plates, 100pcs/Box	100

PCR 384 -well Plate

Picture	Cat#	Description	Packaging (pcs/pk)
	26032020	40uL 384 Well PCR Plate, Full-Skirted, Transparent Plate and White Tube, Black Printing, Cut Corner A24, 10pcs*5bags/Carton	50
	26032021	40uL 384 Well PCR Plate, Full-Skirted, Transparent Plate and Transparent Tube, Black Printing, Cut Corner A24, P24, 10pcs*5bags/Carton	50
	26032022	40uL 384 Well PCR Plate, Full-Skirted, White Plate and Transparent Tube, Cut Corner A24, P24, 10pcs*5bags/Carton	50
	26032023	40uL 384 Well PCR Plate, Full-Skirted, White Plate and White Tube, Cut Corner A24, P24, 10pcs*5bags/Carton	50
	26032024	40uL 384 Well PCR Plate, Full-Skirted, Black Plate and White Tube, Cut Corner A24, 10pcs*5bags/Carton	50